


REF.	10.MC.7101
Descrição / Description Description / Descripción	Mesa de Recém Nascidos, (Berço em Inox) Desmontável / New Born Table, (Stainless Steel Cradle) Removable / Tableau Nouveau-nés (berceau Inox) amovible Mesa de Recién nacidos, (Cuna en acero inoxidable) Desmontable
Material Artigo / Material Article Materiél Article / Material artículo	Aço Inoxidável + Acrílico / Stainless Steel + Acrylic Acier inoxydable + acrylique / Acero Inoxidable + Acrílico
Dimensões Artigo / Dimensions Article Article Dimensions / Dimensiones Artículo	832x455x990mm
Peso máximo paciente / Maximum patient weight / Poids maximum patient / Peso máximo paciente	10Kg
Peso Líquido / Net Weight / Poids Net / Peso Neto	19kg
Peso Bruto / Gross Weight / Poids Brut / Peso Bruto	20kg
Volume Bruto / Gross Volume / Volume Brut / Volumen Bruto	0,539m ³



REF.	10.MC.7102
Descrição / Description Description / Descripción	Berço Pediátrico inclinado, pintado Dim: 680x380x950mm / Pediatric Cradle inclined, painted Dim: 680x380x950mm / Berceau pédiatrique incliné, peint Dim: 680x380x950mm / Cuna Pediátrica inclinada, pintado Dim: 680x380x950mm
Material Artigo / Material Article Materiél Article / Material artículo	Estrutura em Aço acabamento em epoxy / Steel structure Epoxy finishing / Structure en acier finition en epoxy / Estructura en Acero acabado en epoxy
Dimensões Artigo / Dimensions Article Article Dimensions / Dimensiones Artículo	780x512x955mm
Peso máximo paciente / Maximum patient weight / Poids maximum patient / Peso máximo paciente	10Kg
Peso Líquido / Net Weight / Poids Net / Peso Neto	15kg
Peso Bruto / Gross Weight / Poids Brut / Peso Bruto	16kg
Volume Bruto / Gross Volume / Volume Brut / Volumen Bruto	0,342m ³