



REF.	10.AD.1412	10.AD.1402
Descrição Description Description Descripción	Degrau Triplo c/ tela anti-derrapante Triple access step with anti slip screen Marchepied 3 marches avec toile antidérapant Banquillo triplo con pantalla anti deslizamiento	
Material Artigo Material Article Materiel Article Material artículo	Aço com acabamento Epoxy / Steel with Epoxy Finishing / Acier avec finition époxy / Acero con acabado Epoxy	Aço Inoxidável Stainless steel Acier inoxydable Acero inoxidable
Dimensões Artigo / Dimensions Article Article Dimensions / Dimensiones Artículo	660x380x610mm	
Peso máximo admissível / Admissible maximum weight / Poids maximum admissible / Peso máximo admisible	120kg	
Peso Líquido / Net Weight Poids Net / Peso Neto	10,2kg	10,0kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	10,5kg	10,3kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,235m <sup>3</sup>	



REF.	10.AD.1430	10.AD.1431
Descrição Description Description Descripción	Degrau Duplo c/ tela anti-derrapante Double access step with anti slip screen Marchepied 2 marches avec toile antidérapant Banquillo doble con pantalla anti deslizamiento	
Material Artigo Material Article Materiel Article Material artículo	Aço com acabamento Epoxy / Steel with Epoxy Finishing / Acier avec finition époxy / Acero con acabado Epoxy	Aço Inoxidável Stainless steel Acier inoxydable Acero inoxidable
Dimensões Artigo / Dimensions Article Article Dimensions / Dimensiones Artículo	500x380x400mm	
Peso máximo admissível / Admissible maximum weight / Poids maximum admissible / Peso máximo admisible	120kg	
Peso Líquido / Net Weight Poids Net / Peso Neto	7,2kg	6,6kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	7,5kg	7kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,087m <sup>3</sup>	



REF.	10.AD.1432	10.AD.1433
Descrição Description Description Descripción	Degrau Simples c/ tela anti-derrapante Simple access step with anti slip screen Marchepied simple avec toile antidérapant Banquillo simple con pantalla anti deslizamiento	
Material Artigo Material Article Materiel Article Material artículo	Aço com acabamento Epoxy / Steel with Epoxy Finishing / Acier avec finition époxy / Acero con acabado Epoxy	Aço Inoxidável Stainless steel Acier inoxydable Acero inoxidable
Dimensões Artigo / Dimensions Article Article Dimensions / Dimensiones Artículo	485x380x265mm	
Peso máximo admissível / Admissible maximum weight / Poids maximum admissible / Peso máximo admisible	120kg	
Peso Líquido / Net Weight Poids Net / Peso Neto	3,8kg	3,6kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	4,0kg	3,8kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,051m <sup>3</sup>	