


REF.	10.CD.1305
Descrição / Description Description / Descripción	Cadeira Banho Fixa inox - Altura regulável Stainless Steel fixed shower seat - Height adjustable Siège de douche fixe en acier inoxydable - Hauteur réglable / Silla de baño Fija inox - Altura ajustable
Material Artigo / Material Article Materiel Article / Material articulo	Aço Inoxidável / Stainless steel Acier inoxydable / Acero inoxidable
Dimensões Assento / Seat Dimensions Dimensions du siège / Dimensiones asiento	460x405x(346-931)mm
Dimensões Artigo / Dimensions Article Article Dimensions / Dimensiones Artículo	920x750x1050mm
Peso máximo paciente / Maximum patient weight / Poids maximal patient Peso máximo usuario	135kg
Peso Líquido / Net Weight Poids Net / Peso Neto	26kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	27kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,815m <sup>3</sup>



REF.	10.CD.1300
Descrição / Description Description / Descripción	Cadeira de Banho Móvel Inox - Altura Regulável em 3 posições / Stainless Steel mobile shower Chair - Adjustable Height in 3 positions / Siège de douche mobile en acier inoxydable - Hauteur réglable en 3 positions / Silla de baño móvil Inox - Altura regulable en 3 posiciones
Material Artigo / Material Article Materiel Article / Material articulo	Aço Inoxidável / Stainless steel Acier inoxydable / Acero inoxidable
Dimensões Assento / Seat Dimensions Dimensions du siège / Dimensiones asiento	460x405x(491/571/651)mm
Dimensões Artigo / Dimensions Article Article Dimensions / Dimensiones Artículo	600x600x900mm
Peso máximo paciente / Maximum patient weight / Poids maximal patient Peso máximo usuario	135kg
Peso Líquido / Net Weight Poids Net / Peso Neto	18kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	19kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,583m <sup>3</sup>